

TSG 23cu.ft./650L Laboratory Refrigerator, 5C Setpoint, Glass Door, 208-230V 50/60Hz, NEMA 6-15

Thermo Scientific TSG Series Glass Door Laboratory Refrigerators offer quiet, intuitive operation with excellent temperature performance and an outstanding warranty, meeting the needs of a wide variety of professional users who need consistent, reliable performance without sacrificing usability or sustainability. All TSG Series Laboratory Refrigerators utilize SNAP and F-Gas compliant hydrocarbon refrigerants and are configured to provide excellent temperature performance while minimizing power consumption, all delivered from our zero waste to landfill facility in Asheville, North Carolina.

The drive for sample integrity and security

- Designed to safely store media, reagents, and other laboratory-grade products requiring a temperature range of 3°C to 7°C
- Easy to use microprocessor control with interior light, setpoint security and alarms
- Standard door locks provide comfort and assurance



ThermoFisher
SCIENTIFIC