

**Solaris 4000 R Large Incubated and Refrigerated Benchtop
Orbital Shaker with 18"x18" platform with High Adhesion Mat
PN 75004127 280x170mm and Low Adhesion Mat PN 75004111
200x200mm**

Thermo Scientific Solaris Incubated and Refrigerated Benchtop Shakers combine an ergonomic design with key feature enhancements like Peltier chillers, a sealed electronic system, programmable display, and an easy to clean platform design. The triple eccentric drive mechanism, engineered to perform with a reduced noise level, is built with proven technology to support heavy loads, always maintaining consistent performance. The belt is maintenance-free and engineered to last the life of the shaker, and the shaker mechanism itself comes with an industry leading warranty.

