

Multidrop Combi SMART+

Meet your lab's specific needs with exceptional flexibility and performance with the easy-to-use Thermo Scientific™ Multidrop™ Combi Reagent Dispenser, designed to meet the varied requirements of pharmaceutical and biotech laboratories. It combines the most versatile reagent dispensing functions with excellent performance for drug discovery or genomic and proteomic assays.

The Multidrop Combi Reagent Dispenser covers a wide variety of plates and high volume ranges, while offering fast dispensing and intensive operation.

- It facilitates precise dispensing in the range of 0.5 to 2500 μL , ensuring reproducible test data.
- Supports microplates from 6 to 1536 wells, with plate heights from 5 to 50 mm.
- Icon-based graphical viewer that simplifies use and programming, even without training.
- Minimal dead volume and recoil levels reduce reagent cost
- Easy-to-use Thermo Scientific FILLit software increases flexibility and functionality
- Full robotic compatibility improves workability
- It uses autoclavable and removable eight-channel dispenser cassettes that are included as standard across the Multidrop range.

