

**Solaris 4000 I Large Incubated Benchtop Orbital Shaker with
18"x18" platform with Flask Clamp Starter Kit PN 30129 (1)
Microplate Clamp and (2) 125ml, (4) 250ml, (4) 500ml, (2) 1L,
and (2) 2L Flask Clamps**

Thermo Scientific Solaris Incubated and Refrigerated Benchtop Shakers combine an ergonomic design with key feature enhancements like Peltier chillers, a sealed electronic system, programmable display, and an easy to clean platform design.

The triple eccentric drive mechanism, engineered to perform with a reduced noise level, is built with proven technology to support heavy loads, always maintaining consistent performance. The belt is maintenance-free and engineered to last the life of the shaker, and the shaker mechanism itself comes with an industry leading warranty.

- Create, download, and share up to 99 programs across multiple units for better efficiency and reproducibility in the lab
- Full color, touch screen display is easy-to-read, allowing for clear view of the speed, temperature, and run times from a distance

