

WiFi ULT/Walk-In Freezer

Thermo Scientific™ InSight is a wireless monitoring solution that has been designed to proactively monitor your cold storage equipment. InSight consists of wireless modules (motes), sensors and software tools that allow you to monitor your cold storage equipment's temperature, energy consumption, door access, reliability, predictive failure analytics, as well as other parameters using Wi-Fi or proprietary wireless protocols. Alarms are issued automatically to users' email and mobile devices when setpoint breaches are detected in the system.

InSight Analytics:

InSight Analytics, powered by KLATU™, is a system designed to collect, analyze and make efficiency and operating cost recommendations from multiple asset inputs and events such as:

- Cabinet temperature
- Electrical use and energy consumption
- Door openings
- Ambient temperature

