

## Orion Star A323 dissolved oxygen portable meter with batteries

Get high accuracy and premium performance where you need it with the Thermo Scientific™ Orion Star™ A323 Portable Dissolved Oxygen and RDO Meter. This meter has been designed for a wide range of field applications and dissolved oxygen and temperature testing. Evaluate samples in the most challenging locations with this waterproof meter with IP67-rated housing. Use a polarographic DO probe or optical resistant dissolved oxygen (RDO) and calibrate with water-saturated air, air-saturated water, custom value (Winkler), or zero-point. Record up to 5000 sets of time- and date-stamped data points and easily transfer the records to a printer or computer.

Take measurements at the point of origin with our line of superior quality portable meters. The Thermo Scientific Orion Star A323 Portable RDO/DO Meter offers an information-rich display and easy operation. This rugged, waterproof portable meter is ready to take on the elements with you. The Orion Star A323 Portable RDO/DO Meters are reliable and easy-to-use field meters that offer advanced functions and features, including 5000-point date stamping with user ID and sample ID options, high/low limit alarms and calibration and password-protected measurement methods.

- 



**ThermoFisher**  
SCIENTIFIC