

Orion Star A321 pH portable meter with batteries

Take measurements at the point of origin with our line of superior quality portable meters. The Thermo Scientific Orion Star A321 Portable pH Meter offers an information-rich display and easy operation.

This rugged, waterproof portable meter is ready to take on the elements with you. The Orion Star A321 portable pH meters are reliable and easy-to-use field meters that offer advanced functions and features, including 5000-point date stamping with user ID and sample ID options, high/low limit alarms and password-protected calibration and measurement methods.

- Measure relative pH, mV, mV, or ORP with temperature
- Perform pH calibrations of up to five points with automatic buffer recognition of NIST and DIN buffer groups with the option to manually enter custom buffer values
- Calibration edit option allows individual points to be fixed without the need for a complete recalibration
- Informative, easy-to-read graphical display with backlight and clear on-screen instructions, simplifying operation

