

STAR T940 ALL-IN-ONE TITRATOR

Use the Orion Star T940 All-in-One Titrator for flexible pH, redox, and ion concentration titrations, including equivalence point titrations, preset pH or mV endpoint titrations, and Multiple Known Addition (MKA) mode for automatic, known addition of multiple ions. Combine the titrator option with any of our high-end ROSS pH or redox and Orion ion selective electrodes or select one of our pre-configured titrator and electrode kits for a bundle discount.

- Accurate and reproducible results: Automating the titration workflow allows analysis to be performed with the same configuration parameters and calculations every time
- Accelerated setup time: Set up and save electrode, titrator, and titration parameters once with the help of methods
 - Better lab performance: Start titration and let the instrument run unattended
- Reduced costs and downtime: Easily replace titrator parts to minimize downtime for maintenance; Plus, durable electrodes reduce cost of ownership

