

TSB Ref 400 Glass 230V/50Hz UK

Thermo Scientific TSB Series Blood Bank Refrigerators are designed to meet stringent requirements for the storage of blood and blood-related products in accordance with Directive 93/42/EC - Class IIa.

TSB series blood bank coolers are specifically designed to meet the strict standards of Directive 93/42/EC - Class IIa on medical devices for blood storage and related products.

- Wide variety of sizes, from 140 l undercounter to 1500 l double door models.
- Digital controller with touch screen and integrated data logging.
- Very quiet operation <50 dBA.
- Triple leaf glass doors.
- Interior and exterior painted in white.
- Temperature alarms, door open and power cut (with backup battery).
- Low consumption with real-time energy consumption monitoring.
- Low heat emission.
- Rear access port.
- Electronic lock with password protection.



ThermoFisher
SCIENTIFIC