

Precision 19L Circulating Bath

Thermo Scientific Precision™ Water Circulation Baths™ are perfect for cases where uniformity and temperature control are especially critical, such as when working with enzymes and in serological applications. Available in three different models, these high-performance baths feature microprocessor control and include a removable stainless steel two-slope cover, a diffuser tray and a rubber ducky.

- Achieve greater temperature uniformity with perimeter-directed water flow
- Easy cleaning and maintenance with the internal coil-free design
- Optimize scheduling with automatic on/off timers
- Enjoy a higher space to store workshop labware with the new hinged lid
- Help prevent bath overheating and damage with protection in the event of low fluid level
- Easy operation and monitoring with icon-based graphical display
- Protect your work with audible alarms

