

Precision 15L Dubnoff SWB

Thermo Scientific™ Precision Stirred Water Baths™ support a wide range of critical life sciences and quality assurance and control applications—from heating fragile reagents to tissue culture and genetic sequencing. Easy cleaning and maintenance with the new coil-free interior. The Precision Shallow Stirred Water Bath has a removable tray that allows it to be used with smaller sample containers. Precision Dubnoff Stirred Water Baths are specifically designed for applications that require samples to be incubated in a controlled environment.

- Easy cleaning and maintenance with the coil-free design
- Optimize scheduling with automatic on/off timers
- Help prevent bath overheating and damage with protection in the event of low fluid level
- Easily save frequently used parameters with four temperature and stirring speed presets
- Help protect your work with audible alarms
- Easy operation and monitoring with icon-based graphical display
- Enjoy a higher space to store workshop labware with the new hinged lid
- Includes an innovative stirrer motor

