

## Precision 19L Coliform

Make fecal coliform testing easier with Thermo Scientific™ Precision™ coliform water baths, designed specifically for this application. The baths are pre-set to optimal temperatures of 35.0, 41.5, 44.5 and 45.5 °C for ease of use. All baths are calibrated at the factory, with subsequent field calibration from the front panel without the need for tools. Benefit from outstanding corrosion and chemical resistance thanks to the epoxy-coated outer surface and easily clean the chamber with stainless steel interior.

- Temperature controller provides uniformity of  $\pm 0.05^{\circ}\text{C}$  and stability of  $\pm 0.1^{\circ}\text{C}$  with stainless steel two-part cover
- Achieve excellent temperature uniformity with perimeter-directed water flow
- Easy cleaning and maintenance with the coil-free design
- Optimize scheduling with automatic on/off timers
- Help prevent overheating of the bath with protection in case of low liquid level
- Easy operation and monitoring with icon-based graphical display
- Enjoy a higher space to store workshop labware with the hinged lid
- Protect your work with audible alarms

