

TSX REF BLOOD 50 cf 230v/50hz

Thermo Scientific TSX High-Performance Blood Bank™ Refrigerators are designed with features to meet sample protection and sustainability goals for whole blood and blood component storage. The series' V-drive technology has been designed to provide uniform temperature that continuously adapts to usage patterns and offers significant energy savings without compromising protection.

- Designed to meet stringent AABB and FDA requirements for the storage of whole blood and blood components
- Stores up to 385 blood bags (450 mL whole blood)
- It is a quiet device, with a sound emission level of 52 dBA, which allows the end user to work next to the refrigerator without any distractions
- The unique V-drive technology is designed to detect usage patterns, such as door openings, if higher compressor speed is needed, as well as periods of stability when the compressor runs at a lower speed; All this allows energy savings without jeopardizing protection

