

CORROSION-RESISTANT VACUUM PU

Achieve maximum solute recovery with Thermo Scientific Deep Vacuum™ Oil Pumps for Vacuum Concentrators, designed to work with all Thermo Scientific™ SpeedVac™ concentration and freeze drying systems.

- Three-positional gas flow control optimizes pump performance with all applications
- Large side glass sight glass makes it easy to check for proper oil levels
- All rotary vane oil pumps come standard:
 - Muffler (OME190 Oil Vapor Eliminator) Two Oil Filters (MF 190) Quick Drain Kit

