

TSX FRZ MANUAL 23cf 208v/60hz

- Designed for the safe storage of medical and laboratory products, including vaccines, pharmaceuticals, media, and reagents that need to be at a temperature between 2°C and 8°C
- It is a quiet device, with a sound emission level of 52 dBA, which allows the end user to work next to the refrigerator without any distractions
- The unique V-drive technology is designed to detect usage patterns, such as door openings, if higher compressor speed is needed, as well as periods of stability when the compressor runs at a lower speed; All this allows energy savings without jeopardizing protection
- With up to 40% lower energy consumption than conventional refrigerant models*
- Natural R290 refrigerants that enable TSX Series devices to meet the U.S. Environmental Protection Agency's Significant New Alternatives Policy (SNAP) standard
- Manufactured in an award-winning waste-to-landfill facility

