

## TSX REF PHARM 23cf 208v/60hz

- Thermo Scientific™ TSX Series high-performance pharmacy refrigerators are designed with features to meet sample protection and sustainability goals for the storage of pharmaceuticals, vaccines, chemotherapy, and other types of pharmaceuticals and medical products that need to be between 2° and 8°C. The series' V-drive technology is designed to provide a uniform temperature that continuously adapts to usage patterns and offers significant energy savings without compromising protection.
- Designed for the safe storage of pharmaceuticals, vaccines, chemotherapy and other medical products that need to be between 3°C and 7°C
- It is a quiet device, with a sound emission level of 52 dBA, which allows the end user to work next to the refrigerator without any distractions
- The unique V-drive technology is designed to detect usage patterns, such as door openings, if higher compressor speed is needed, as well as periods of stability when the compressor runs at a lower speed; All this allows energy savings without jeopardizing protection
- With up to 40% lower energy consumption than conventional refrigerant models\*
- Natural R290 refrigerants that enable TSX Series devices to meet the U.S. Environmental Protection Agency's Significant New Alternatives Policy (SNAP) standard

