

TSX FRZ ENZYME 23cf 230v/50hz

- Thermo Scientific™ TSX Series High-Performance -20°C Enzyme Freezers with Manual Defrost are designed with features to meet sample protection and sustainability goals for enzyme storage that needs to be at -20°C. The V-drive technology in this series is designed to provide a uniform temperature that continuously adapts to the user's patterns and enables significant energy savings without compromising protection.
- It is a quiet device, with a sound emission level of 52 dBA, which allows the end user to work next to the refrigerator without any distractions
- The unique V-drive technology is designed to detect usage patterns, such as door openings, if higher compressor speed is needed, as well as periods of stability when the compressor runs at a lower speed; All this allows energy savings without jeopardizing protection
- With up to 50% lower energy consumption than conventional refrigerant models*
- Natural R290 refrigerants that enable TSX Series devices to meet the U.S. Environmental Protection Agency's Significant New Alternatives Policy (SNAP) standard
- Manufactured in an award-winning waste-to-landfill facility
- Temperature range from -35° to -15°C, factory default to -20°C



ThermoFisher
SCIENTIFIC