

TSX REF BLOOD 12cf 208v/60Hz

- Thermo Scientific TSX High-Performance Blood Bank™ Refrigerators are designed with features to meet sample protection and sustainability goals for whole blood and blood component storage. The series' V-drive technology has been designed to provide uniform temperature that continuously adapts to usage patterns and offers significant energy savings without compromising protection.
- Designed to meet stringent AABB and FDA requirements for the storage of whole blood and blood components
- Stores up to 385 blood bags (450 mL whole blood)
- It is a quiet device, with a sound emission level of 52 dBA, which allows the end user to work next to the refrigerator without any distractions
- The unique V-drive technology is designed to detect usage patterns, such as door openings, if higher compressor speed is needed, as well as periods of stability when the compressor runs at a lower speed; All this allows energy savings without jeopardizing protection
- With up to 40% lower energy consumption than conventional refrigerant models*
- Natural R290 refrigerants that enable TSX Series devices to meet the U.S. Environmental Protection Agency's Significant New Alternatives Policy (SNAP) standard



thermoFisher
SCIENTIFIC