

TSX FRZ LAB 12cf 230v/50Hz

- Thermo Scientific™ TSX Series High-Performance -30°C Cool-Defrost Freezers are designed with features to meet sample protection and sustainability goals for -30°C storage of clinical and laboratory products. The series' V-drive technology is designed to provide a uniform temperature that continuously adapts to usage patterns and offers significant energy savings without compromising protection.
- It is a quiet device, with a sound emission level of 52 dBA, which allows the end user to work next to the refrigerator without any distractions
- The unique V-drive technology is designed to detect usage patterns, such as door openings, if higher compressor speed is needed, as well as periods of stability when the compressor runs at a lower speed; All this allows energy savings without jeopardizing protection
- With up to 50% lower energy consumption than conventional refrigerant models*
- Natural R290 refrigerants that enable TSX Series devices to meet the U.S. Environmental Protection Agency's Significant New Alternatives Policy (SNAP) standard

