

**Solaris 4000 R Large Incubated and Refrigerated Benchtop
Orbital Shaker with 18"x18" platform with Flask Clamp Starter
Kit PN 30129 (1) Microplate Clamp and (2) 125ml, (4) 250ml, (4)
500ml, (2) 1L, and (2) 2L Flask Clamps**

- Create, download, and share up to 99 programs across multiple units for better efficiency and reproducibility in the lab
- Full color, touch screen display is easy-to-read, allowing for clear view of the speed, temperature, and run times from a distance
- User interface has 10 optional languages
- Built-in orbital calculator aids in converting protocols designed with different orbits and estimating the set speed needed to obtain similar results
- Efficient and Reliable Operation
- Reliable, heavy-duty triple eccentric drive mechanism provides timed or continuous 24-hour operation to allow smooth and uniform shaking of large or uneven loads
- Unbalanced load sensor stops platform motion when excess vibration is detected and visual/audible alarms signal until condition is corrected

