

TSG 23cu.ft./650L Laboratory Refrigerator, 5C Setpoint, Glass Door, 115V/60Hz, NEMA 5-15P

- Thermo Scientific TSG Series Glass Door Laboratory Refrigerators offer quiet, intuitive operation with excellent temperature performance and an outstanding warranty, meeting the needs of a wide variety of professional users who need consistent, reliable performance without sacrificing usability or sustainability. All TSG Series Laboratory Refrigerators utilize SNAP and F-Gas compliant hydrocarbon refrigerants and are configured to provide excellent temperature performance while minimizing power consumption, all delivered from our zero waste to landfill facility in Asheville, North Carolina.
- The drive for sample integrity and security
- Designed to safely store media, reagents, and other laboratory-grade products requiring a temperature range of 3°C to 7°C
- Easy to use microprocessor control with interior light, setpoint security and alarms
- Standard door locks provide comfort and assurance



ThermoFisher
SCIENTIFIC