

TSG 23cu.ft./650L Laboratory Refrigerator, 5C Setpoint, Solid Door, 208-230V 50/60Hz, NEMA 6-15

- Designed to safely store pharmaceuticals & vaccines requiring a temperature range of 2°C to 8°C
- Easy to use microprocessor control with interior light, setpoint security and alarms
- Standard door locks provide comfort and assurance
- Enabling compliance
- Tested and Certified to NSF/ANSI 456 Vaccine Storage standard and supports the US CDC, WHO, and other global vaccine storage guidelines and recommendations
- Tested and certified to conform with the design and safety standards of UL 61010, and CE Directive 2011/65/EU and 2015/863

