

TSG 30cu.ft./850L NSF 456 Pharmacy Refrigerator, 5C Setpoint, Glass Door, 208-230V 50/60Hz, NEMA 6-15

- Designed to safely store media, reagents, and other laboratory-grade products requiring a temperature range of 3°C to 7°C
- Easy to use microprocessor control with interior light, setpoint security and alarms
- Standard door locks provide comfort and assurance
- Enabling compliance
- Tested and certified to conform with the design and safety standards of UL 61010, and CE Directive 2011/65/EU and 2015/863
- Designed for an ergonomic and comfortable experience
- Ergonomically designed door handle
- Convenient, easy to clean and adjustable epoxy coated full and half shelves
- Whisper quiet with a sound level of 52 dBA, enabling end users to work alongside the refrigerator without distraction
- Bright, LED interior lighting
- Self-closing multipaned glass doors with 90° stop to assist with inventory loads

