

EL CAST BRNZ 1"NOM CXM

- Confidently perform silica analysis with the Thermo Scientific™ Orion™ 8030cX Silica Analyzer. Designed to provide accurate, effortless water analysis, the Orion 8030cX specializes in online measurement of reactive silica. Powered by advanced colorimetric technology, the Orion 8030cX delivers premium precision, greater productivity and enhanced durability. Supporting ultrapure and process water applications, the Orion 8030cX conveniently delivers precise, reliable measurement you can depend on. Unleash the possibilities and conquer the freedom of fast, confident decision making.
- Confident Results
- Low Detection Limit (0.5 ppb)
- High Accuracy (≤ 1 ppb or 1% for range 0-500 ppb)
- Patented Fluid Delivery & Quantification System
- Reduced Phosphate Interferences
- Low Cost Ownership
- 100 Days of Continuous Run Time (w/ 2L reagent and 15-minute measurement frequency)
- Low Reagent Consumption (200 μ L per measurement)
- Separate Chemical Waste Collection (decreased disposal expenses)



ThermoFisher
SCIENTIFIC