

STAR T900 SERIES DISPENSE KIT

- Get accessories and replacement parts for the Orion Star T910 pH Titrator, Orion Star T920 Redox Titrator, Orion Star T930 Ion Titrator, and Orion Star T940 All-in-One Titrator. Designed to increase functionality and reduce maintenance downtime.
- Choose from different interchangeable burette sizes (10 mL, 20 mL, or 50 mL) to customize your titrator
- Bottle holder helps secure the titrator to prevent spillage and allows for flexible installation
- Color-coded tubes and connectors simplify tube installation and maintenance

