

CONV MID FRZ 01 W TS

- Manual Defrost countertop freezer allows for easy transportation and placement within any laboratory or clinical setting, with a countertop footprint about the size of a microwave. These freezers increase user efficiency and lessen physical strain, as there is no bending or squatting required. Top-mounted controls, reversible doors and an ergonomic handle contribute to its modern design, which allows for a convenient, practical approach to a busy day. Multiple GPFs in a medical facility allow for the convenient storage of daily doses of vaccines and pharmaceuticals requiring frozen storage in multiple locations.
- Designed to safely store medical and laboratory storage, including vaccines and pharmaceuticals requiring a temperature of -20°C
- Manual defrost
- External display provides visual confirmation of the internal temperature
- Audible and visual hi/low temperature alarm will alert personnel if temperature levels rise or fall to an undesirable degree
- Compatible with wireless monitoring systems and remote alarm contacts
- Standard door locks provide comfort and assurance
- Insulation and natural refrigerants meet the standards established by the U.S. Environmental Protection Agency's Significant New Alternatives Policy (SNAP)
- Quiet performance with less than 40 decibels of sound
- Ergonomic design including reversible doors and comfortable door handle
- cULus

