

P82DS ADAPTER PLATE WITH

- Ensure smooth, flat belts with crisp resolution with Thermo Scientific™ Owl™ P82 Dual Gel System replacement parts. For use with the Owl™ P82 Dual Gel System; Parts include replacement plates, gaskets, and shims.
- Purchase the necessary spare parts to ensure that flat results are achieved even with irregular or scratched resolutions.

For use with (equipment)P82 Vertical Electrophoresis System

Product LineOwl

TypeAdapter Plate

Unit SizeEach

