

PLATFORM SUSPENSION ARMS

Use Thermo Scientific™ accessories for MaxQ™ 7000 Water Bath Stirrers to increase the utility of MaxQ7000 Water Bath Stirrers. Exceptional innovation, reliability, and accuracy make Thermo Scientific MaxQ7000 shakers, platforms, and accessories the optimal job choice for any laboratory.

For use with (application)for near-ambient temperature and below (requires a circulating chiller)

For use with (equipment)MaxQ 7000

Product TypeCooling coil

Unit SizeEach

