

SUPER NUOVA S 4POS 220-240V

- Perform solubility, digestion, titration, and enzymatic studies with Thermo Scientific Super-Nuova™ Multi-Position™ Digital Shakers. With four individually controlled stirring positions, ceramic tops and 4-function memory, these stirrers are designed to offer a more robust coupling and stirring control with advanced safety features. Ideal for multi-user labs.
- Key features
- The Thermo Scientific™ Stir Trac™ Stirring Control offers improved low-speed stirring, smooth speed, and powerful magnetic coupling, reducing the risk of stir bar decoupling
- Microprocessor control ensures precise and stable stirring speed adjustments
- Stir Trac braking feature completely stops the stirring bar for immediate flask removal
- The single control knob controls the speed; Lock function prevents accidental changes
- Digital rpm display
- Cleans easily and resists alkalis and acids
- One-piece white reflective ceramic surface promotes visibility of samples
- Adjustable stirring in 1 rpm units
- RS-232 port allows you to transmit time counting and stirring speed (adjustment) data directly to your PC

