

STAR T910 PH TITRATOR WITHOUT ELECTRODE

- Use the Orion Star T910 pH Titrator for specific titrations of acids and bases, including titratable acidity of juices and wines, acidity of food products, alkalinity of waters, alkalinity and acidity of consumer products, total acid number (TAN), and total base number (TBN). Titration techniques include equivalence point titrations and predefined pH endpoint titrations to enable versatile sample analysis. Combine the titrator option with any of our high-end ROSS pH electrodes or select one of our pre-configured titrator and electrode kits for a bundled discount.
- •Accurate and reproducible results: Automating the titration workflow allows analysis to be performed with the same configuration parameters and calculations every time
- •Accelerated setup time: Set up and save electrode, titrator, and titration parameters once with the help of methods
- •Better lab performance: Start titration and let the instrument run unattended
- •Reduced costs and downtime: Easily replace titrator parts to minimize downtime for maintenance; Plus, durable electrodes reduce cost of ownership
- •Securely document test results: View live results in the summary table and 100-point

