

FSE PH METER SET; 3-IN-1 ELECT

- Designed for convenient routine pH measurement in the lab: user-friendly, stores data sets in non-volatile memory and meter-attachable electrode stand. Efficient sample reading and recording with flexible display options: display both pH and temperature or mV and temperature simultaneously.
- The essential functions needed for routine measurements
- One to five point pH calibration
- Automatic buffer recognition with selection of two pH buffer sets
- pH slope and offset display
- One point mV calibration
- One point ATC temperature probe calibration
- Automatic or manual temperature compensation
- Ready icon to indicate measurement stability
- Electrode condition icon to indicate status
- 100 data set storage with non-volatile memory

