

TSG Ref 400 Glass 230V/50Hz UK

- Thermo Scientific TSG Series General Purpose™ Laboratory Refrigerators are designed with your needs in mind. The TSG series is ideal for the demanding daily use requirements of research, clinical, and industrial laboratories that store reagents, culture media, proteins, and samples used daily in the laboratory.
- TSG Series General Purpose Refrigerators are designed with your laboratory in mind, combining the performance, security, and quality required for everyday laboratory applications.
- Wide range of sizes from 400 L to 1500 L double door
- Glass door (triple pane) and solid door configurations
- Whisper quiet operation < 50 dBA
- White painted interior and exterior
- Adaptive Defrost
- Temperature, door ajar, and power failure (with battery backup) alarms.
- Low heat output
- Rear access port
- Keyed Lock
- LED interior lighting
- Green hydrocarbon refrigerants
- Remote alarm contacts



ThermoFisher
SCIENTIFIC