

TSG Ref 700 Glass 230V/50Hz UK

- Thermo Scientific TSG Series General Purpose™ Laboratory Freezers are designed with your needs in mind. The TSG series is ideal for the demanding daily use requirements of research, clinical, and industrial laboratories that store reagents, culture media, proteins, and samples used daily in the laboratory.
- TSG Series Freezers are designed with your laboratory in mind, combining the performance, security, and quality required for everyday laboratory applications.
- Wide range of sizes from solid door 400 L to 1500 L double door.
- Whisper quiet operation < 50 dBA
- White painted interior and exterior
- Adaptive Defrost
- Temperature, door ajar, and power failure (with battery backup) alarms.
- Low heat output
- Rear access port
- Keyed Lock
- LED interior lighting
- Green hydrocarbon refrigerants
- Remote alarm contacts



ThermoFisher
SCIENTIFIC