

Precision 5 & 10L Dual GP Bath

- Support a wide range of applications with Thermo Scientific Precision General Purpose™ Water Baths™, designed to maintain water temperature from room temperature to 100°C. These solid, high-performance baths range from 2 to 28 L, including shallow models. The over-temperature safety circuit has been designed to prevent thermal run-out, while new automatic on/off timers allow for optimised operating hours. Benefit from outstanding corrosion and chemical resistance thanks to the epoxy-coated outer surface and easily clean the tank with stainless steel interior.
- Reclaim valuable bench space thanks to its small footprint compared to previous models
- Advanced microprocessor controller designed for greater functionality
- Protect your work with audible alarms
- Easily save the most commonly used parameters with four temperature presets
- Easy operation and monitoring with icon-based graphical display
- Improved hinged design allows for open lid work and easier lid removal
- Help prevent bath damage and overheating with protection in the event of low liquid level
- Simplify operation with quick-disconnect drain on 10, 20, 28L and dual models

