

ULT; PKG; TDE40040LA+Racks

- Thermo Scientific™ TDE -40°C ultra-low temperature freezer packages include one ULT freezer and two shelves of sliding drawer racks and boxes.
 - Package includes:
 - One Thermo Scientific TDE -40°C Ultra-Low temperature freezer
 - 2 shelves of sliding drawer racks with 2-inch cryoboxes and 100-count cell dividers. (see product specifications for details)
- Dimensiones exteriores (L × An × Al) 38.5 x 28.3 x 78 in. (97.7 x 71.9 x 198.1cm)
Requisitos eléctricos 115 V, 60 Hz
Admite Cajas (400)
Incluye 400 box capacity freezer, 8 sliding drawer racks, 200 2-inch boxes with 100-count cell dividers
Aislamiento Water-Blown Polyurethane Foam
Interior Acero inoxidable



ThermoFisher
SCIENTIFIC